



European
Funds
Smart Growth



Republic
of Poland

European Union
European Regional
Development Fund



KUKINI™

Special needs stroller

KUKINI™ special needs stroller is lightweight, making it ideal for long walks and running. Its aluminum frame with reflective elements ensures safety, while the rear shock-absorbing system and 16-inch inflatable wheels provide a smooth ride on uneven terrain. A 5-point safety belt with a magnetic fastener secures the child, while the Octagon™ headrest offers comfortable support. The adjustable backrest allows for an optimal sitting position. A quick-release wheel system and compact folding design make storage and transport easy.



Light-
weight

24⁺

24-month
warranty



Stable
construction



Smooth and stepless
adjustment

Colors

Tapicerka HYDRO™

Waterproof and durable



Cleverly designed, lightweight aluminium frame guarantees safety and comfort.



A special bicycle tow bar can turn the tricycle into a bike trailer.



100 - 165 cm

Indicative user
height



40 - 65 kg

Max. user
weight



For more details about
the product, please visit
our website.



This is a medical device.
Use it in accordance with the user manual or the label.

INDICATIONS

Cerebral palsy (CP) – various types and forms | Brain and spinal cord injuries of various etiologies with limb paresis, severe balance disturbances, and coordination impairments | Demyelinating diseases (multiple sclerosis – MS, amyotrophic lateral sclerosis – ALS, and others) | Meningocele | Muscle diseases with paresis and paralysis (including dystrophy, myasthenia gravis, polymyositis) | Genetic syndromes with limb paresis | Degenerative and metabolic diseases of the nervous system | Degenerative and rheumatic joint diseases | Other diseases with limb paresis, paralysis, and musculoskeletal impairments | After injuries, following traumatic brain injuries or spinal cord injuries with limb paresis and paralysis, particularly with quadriplegia (as part of the rehabilitation process or as permanent orthopedic support)

CONTRAINDICATIONS

Observation: Monitor user during standing. Stop therapy if paleness, sweating, fainting, joint overload, or body deformations occur; consult with a doctor or physiotherapist | **Absolute Contraindication:** Critical condition of the user. | **Relative Contraindications:** Resting hypertension > 200/110 mmHg, pulmonary hypertension, heart rhythm disorders, unstable asthma, unstable metabolic issues, moderate heart valve disease. | Conditions like healing fractures or severe osteoporosis that prevent weight-bearing | Incomplete fracture healing (unless approved by a doctor) | Postoperative joint load restrictions | Pain from dislocation, excessive tissue stretch, or pressure; dislocation without pain is not contraindicated | Significant hip/knee flexion or ankle contractures; ensure no excessive pressure or stretch on affected areas.

Manufacturer: Akces-MED Sp. z o.o. | Poland, 36-002 Jasionka | Jasionka 955B